


ORTHO

REGULATION



SUPPLIES



CRAFTSMANSHIP

ORTHODONTIC
LABORATORY
& SUPPLIES
orthoregulation.eu
+31(0)36 52 95 099

ORTHO

REGULATION

CRAFTSMANSHIP

We deliberately choose to manufacture our products in the Netherlands. A dental brace is a custom-made product that is hand-crafted in our own orthodontic laboratory. Changing **customer demands** and high **levels of quality and service** create the need for a targeted and modified product. For instance, certain adjustment options for braces may be a prerequisite. We focus especially on such demands and in this capacity, we are sure to be able to serve you well at our laboratory. The organizational structure at **Ortho Regulation** is characterized by a simple structure with **short lines of communication**.

For over **twenty years**, we have been manufacturing orthodontic removable and fixed appliances, whereby **craftsmanship and enjoyment** related to creation of our products are our priorities. Of course we use the latest proven technology. Our traditional way of working, supported by our ongoing **innovation** with modern techniques, continues to make our profession the most beautiful thing we can think of. We always maintain a critical eye for detail and with this in mind it is possible for us to guarantee the quality you expect from us as specialists; together with a personal approach that works best for you.

Ortho Regulation is an independent laboratory where all of its handcrafted specialized products are made on location and **not delegated to third parties**.

SUPPLIES

Ortho Regulation is not only an **expert** in orthodontic removable and fixed appliances. You can also rely on us for all other orthodontic products for your practice, such as ortho-alginate, brackets, arches, impression systems, instruments etc. Are you unable to find what you are looking for straight away? This is not a problem; we provide **customization** of all necessary equipment and products. We can also deliver any supplies of various brands.

THE FUTURE

Ortho Regulation continuously focusses on **discovering** new possibilities and **trends** in orthodontics and digitalization. Would you prefer to continue to print conventionally, are you changing over to the use of an intra oral scanner or do you share an intra oral scanner in your practice? Are you going to provide digital STL files to your laboratory? How are we going to enter the orders; are we still using the technique slips? Will the order records be digitally saved or are they already delivered digitally? Will we be using 3D printers etc.? The fact that this is becoming the future is certain, but when will this be cost effective and **efficient**?

COMMON GOAL

You are the specialist, you know what you need. Your knowledge and experience increases our expertise and ensures that we continue to learn; even after twenty years. As far as we are concerned, only one thing really counts: a healthy smile for your patient. This is what we want to achieve together with you. For this, we do our utmost day in and day out. We are proud of our proven **stability, continuity and integrity**.



ORTHODONTIC
LABORATORY
& SUPPLIES
www.orthoregulation.eu
+31(0)36 52 95 099

ORTHO

REGULATION

SENDING METHODS

TRADITIONAL METHOD

Day 1



Patiënts impressions taken with cavex orthotrace alginate method.

Day 1-2



Give us a call or email and received an ups ticket by email for sending the impressions to Ortho Regulation. Between 1600 and 1700PM the driver of ups will pick-up the package out of your practice.

Day 2



Our technicians created plaster models for making orthodontic appliance.

Day 5



The orthodontic appliance finish and will be checked.

Day 5-7




The orthodontic appliance are delivered back to your practice.

*We also work with private drivers ask for various possibilities.

DIGITAL METHOD

Day 1



Patient's mouth is scanned very well, using an intra-oral scanner. STL-Files are send to Ortho Regulation

Day 2



Ortho Regulation begin printed 3d models.

Day 5



The orthodontic appliance finish and will be checked for proper fit on the digital 3d printed model.

Day 5-7



The orthodontic appliance are delivered back to your practice.

*We also work with private drivers ask for various possibilities.

ORTHODONTIC
LABORATORY
& SUPPLIES
orthoregulation.eu
+31(0)36 52 95 099